



Yewtwist for pediatrics: one connection at a time, one patient at a time.

Introducing Yewtwist, a unique tool to safely and efficiently disconnect medical connections.

Prevent IV damage with a simple twist

AN ANSWER FOR EVERY DIFFICULT CONNECTION



Central/peripheral IV access lines



Needleless connectors & IV connections



Feeding tubes (ENFit)



Removing blood tubing



Removing spike from IV bag



Nasogastric (NG) tubes



Drainage systems

Improving the management of medical connections with innovative, safe, and efficient tools to empower healthcare professionals.

DESIGNED BY AN ONCOLOGY NURSE



Increase patient safety and comfort using Yewtwist

- **Easily disconnect:** IV and ENFit connections are difficult to disconnect due to tightness and sticky solutions. Common infusions that contribute to stuck connectors include TPN, blood and blood products, and Propofol – all of which require daily line changes.
- **Loosen tight feeding tubes:** These can often get stuck due to over-tightening and sticky solutions.
- **Protect your IV line from damage:** Reduce unnecessary procedures or treatment delays caused by damaged lines.
- **Avoid delays in care:** Delays in treatment can negatively impact patient health and increase both the length of stay and hospital costs.